


## Spheroid Culture in Matrigel on a Pillar Plate

This standard operating procedure (SOP) describes the step-by-step procedure for manually loading spheroids in Matrigel onto a 36PillarPlate and culturing the spheroids using either a 384DeepWellPlate or a 36PerfusionPlate. Please read the protocol carefully before performing the experiments.

### Materials:

- 36PillarPlate (Bioprinting Laboratories Inc., Cat. no. 36-01-00)
- LoadingPlate (Bioprinting Laboratories Inc., Cat. no. 384-03-00)
- 36PetriLid (Bioprinting Laboratories Inc., Cat. no. 36-03-00)
- 384DeepWellPlate (Bioprinting Laboratories Inc., Cat. no. 384-02-00)
- Ultralow attachment (ULA) 384-well plate (S-BIO, Cat. no. MS-9384UZ; FaCellitate, Cat. no. F224384)
- Growth factor reduced Matrigel (Corning, Cat. no. 354230)
- Deep petri dish, 100 mm x 20 mm (Corning, Cat. no. 70165-102)

### Methods:



The overall protocol of spheroid transfer and attachment to Matrigel on the pillar plate using an ultralow attachment (ULA) 384-well plate.

### Spheroid transfer from ultralow attachment (ULA) 384-well plate to 36PillarPlate and spheroid culture in either 384DeepWellPlate or 36PerfusionPlate

#### Formation of spheroids in an ultralow attachment (ULA) 384-well plate

1. Prepare spheroids in a ULA 384-well plate by seeding 500 - 3,000 cells per well and incubate them for 2 - 4 days.

**Note:** *If spheroids form within 2 days, add 40  $\mu$ L of cell suspension to each 384-well for spheroid formation, then proceed with pillar plate insertion and spheroid transfer without removing the existing cell culture medium. If spheroid formation takes longer than 2 days, requiring a medium change, add 80  $\mu$ L of cell suspension to each 384-well. Before pillar plate insertion and spheroid transfer, carefully remove 40  $\mu$ L of the old medium, leaving 35 - 40  $\mu$ L in each well. This step is*

critical to prevent medium overflow during pillar plate insertion and spheroid transfer.

**Note:** If edge effects are observed, add 80  $\mu\text{L}$  of sterile water to the surrounding wells of the ULA 384-well plate to minimize medium evaporation and promote uniform spheroid formation.

2. Inspect cell spheroids in the ULA 384-well plate under a brightfield microscope prior to spheroid transfer to the pillar plate. The typical size of cell spheroids transferred is 100 - 400  $\mu\text{m}$ .

**Note:** The optimal volume of cell culture medium in each well of the ULA 384-well plate for successful spheroid transfer is 35 - 40  $\mu\text{L}$ . An excess of medium can cause overflow during the pillar plate sandwiching process, while an insufficient volume may result in unsuccessful spheroid transfer due to bubble formation.

### **Preparation of Matrigel stock, cell culture medium, and pillar plate**

3. Thaw Matrigel<sup>®</sup> stock overnight by submerging an unopened bottle in a bucket of ice placed in a 4°C refrigerator. Prepare 500  $\mu\text{L}$  aliquots of Matrigel and store them at - 20°C for future use.

4. Thaw Matrigel<sup>®</sup> aliquots overnight in a 4°C refrigerator prior to spheroid transfer.


**Note:** It is important to thaw Matrigel aliquots in advance in a 4°C refrigerator and maintain Matrigel chilled on ice during use since Matrigel starts to solidify above 10°C. Do not freeze and thaw Matrigel aliquots.

5. For cell culture, dispense 70  $\mu\text{L}$ /well of cell growth medium into a 384DeepWellPlate or 800  $\mu\text{L}$ /fluidic channel of cell growth medium into a 36PerfusionPlate. Cover the plate with an appropriate well plate lid and incubate it in a humidified 5%  $\text{CO}_2$  incubator at 37°C for at least 1 hour prior to use.

**Note:** Prewarming the medium helps minimize temperature shock and reduces air bubble formation during plate assembly and culture. Adding an excessive volume of cell culture medium to the 384DeepWellPlate or 36PerfusionPlate may cause overflow after the pillar plate is sandwiched with the well plate. Avoid wetting the bottom of the pillars with culture medium during this process, as it may result in cross-talk or contamination between wells.

6. Hydrate the surface of the pillar plate by inserting two 36PillarPlates into a 36PetriLid placed on a 100 x 20 mm petri dish containing 5 mL of sterile distilled water. Incubate the assembly in a humidified 5%  $\text{CO}_2$  incubator at 37°C for 30 minutes prior to hydrogel loading (**Fig. 1**).

**Note:** Hydrating the surface of the pillar plate in a humid environment is necessary to increase surface hydrophilicity and minimize air bubble entrapment on the pillars after Matrigel loading. Ensure that the pillars are not immersed in water when transferring the assembly to the  $\text{CO}_2$  incubator, as excess water on the pillars may interfere with uniform hydration.



**Figure 1.** Hydration of the pillar plate surface in a 100 x 20 mm petri dish with 5 mL of sterile, distilled water to minimize air bubble entrapment.

### **Preparation of diluted Matrigel**

7. Gently mix 1.5 mL of warm cell culture medium with 1.5 mL of cold Matrigel to obtain a final concentration of 4 - 6 mg/mL Matrigel (i.e., 50% Matrigel).

**Note:** Use diluted Matrigel immediately after preparation. If not used immediately, keep the diluted Matrigel on ice. Do not refreeze or reuse thawed or diluted Matrigel. For long-term organoid culture, or when severe spheroid detachment from the pillar plate is observed (typically after


approximately 5 - 10 days of cancer spheroid culture due to Matrigel degradation by cell-secreted metalloproteases), use 50% Matrigel supplemented with 0.5% alginate. Alginate is resistant to degradation by metalloproteases and helps maintain hydrogel integrity during extended culture.

8. If bubbles form on the pillar plate after hydrogel loading, centrifuge the 50% Matrigel solution at 1,000 rpm for 2 minutes to remove air bubbles generated during the mixing process.

**Note:** Avoiding micro-bubble entrapment during Matrigel mixing is critical to prevent air bubble formation on the pillars after 1-day culture.


#### **Matrigel loading on single pillar plate using a 1 mL pipette tip (to save Matrigel)**

9. Aspirate 250  $\mu$ L of 50% Matrigel using a 1 mL pipette tip for single 36PillarPlate.  
**Note:** It will require 5  $\mu$ L of Matrigel per pillar (at least 180  $\mu$ L per 36PillarPlate).
10. Separate the 1 mL pipette tip from the pipette gently to prevent Matrigel spillage.
11. Using the index finger, block the back opening of the pipette tip to prevent Matrigel overflow while tapping the pillar surface.
12. Gently tap the 1 mL pipette tip containing 50% Matrigel onto the center of the pillar to load the diluted Matrigel while blocking the large back opening of the tip with the index finger (Fig. 2).  
**Note:** Do not touch the bottom surface of the pillar with the pipette tip to avoid damaging the surface coating. Use the hydrated pillar plates within 5 minutes after removal from the humidified petri dish to prevent complete drying of the surface.



**Figure 2.** Loading Matrigel on the pillar plate using a 1 mL pipette tip with Matrigel.


13. Repeat **Step 12** for all pillars.
14. After loading Matrigel on all pillars, remove excess Matrigel by horizontally sliding a 1 mL pipette tip across the pillar surfaces (**Fig. 3**).



**Figure 3.** Scraping excess Matrigel from the pillars using a 1 mL pipette tip.


#### **Matrigel loading on multiple pillar plates using a LoadingPlate**

15. Place a LoadingPlate on a flat surface. Dispense 2 - 2.5 mL of the diluted Matrigel solution into each small block without introducing big bubbles, and spread the solution evenly using the pipette tip (**Fig. 4**).



**Figure 4.** Dispensing 2 - 2.5 mL of diluted Matrigel per block in the LoadingPlate for rapid loading of Matrigel on the pillar plate.


**Note:** It is critical to maintain a sufficient volume of Matrigel in each block of the LoadingPlate; a minimum volume of 2 mL per block is recommended. Matrigel should be distributed uniformly to ensure complete wetting of all pillars. Improper loading of Matrigel onto the pillars during pillar stamping may result in macro-bubble formation on the pillars after spheroid transfer.



**Note:** Do not leave the diluted Matrigel solution on the LoadingPlate for longer than 5 minutes to avoid premature gelation during the stamping process. Because pillar stamping is performed rapidly, it is generally not necessary to place the LoadingPlate containing diluted Matrigel on ice during this step.

16. Stamp the 36PillarPlate on the LoadingPlate and press gently to load diluted Matrigel evenly on the entire pillar plate. Repeat this Matrigel loading step for another pillar plate (**Fig. 5**).

**Note:** Using 2 - 2.5 mL of the diluted Matrigel solution, it is possible to prepare at least four 36PillarPlates (5  $\mu$ L Matrigel per pillar or 180  $\mu$ L per 36PillarPlate) without introducing macro-bubbles on the pillars. For uniform pillar wetting and robust Matrigel loading, gently wiggle the pillar plate during stamping. Add additional Matrigel solution to the LoadingPlate as needed.



**Figure 5.** Stamping the 36PillarPlate onto the LoadingPlate to load Matrigel onto the pillars.

#### **Complete Matrigel gelation prior to spheroid transfer**

17. For minimizing water evaporation during complete Matrigel gelation, insert two 36PillarPlates loaded with diluted Matrigel into a 36PetriLid placed on a 100 x 20 mm petri dish containing 5 mL of sterile, distilled water (**Fig. 1**).
18. Incubate the assembly in a humidified 5% CO<sub>2</sub> incubator at 37°C for 20 minutes to allow complete gelation of the diluted Matrigel.


**Note:** It is critical to minimize water evaporation during Matrigel gelation to ensure proper spheroid transfer. Ensure that the pillars are not immersed in water when transferring the assembly to the CO<sub>2</sub> incubator, as excess water on the pillars may interfere with proper gelation.

### Spheroid transfer and attachment on the pillar plate

- Align one corner pillar of the 36PillarPlate with the corresponding corner well of the ULA 384-well plate. Carefully sandwich the pillar plate with Matrigel onto the ULA 384-well plate containing spheroids. Cover the assembled plates with a 384-well plate lid and quickly invert the sandwiched plates so that the pillar plate faces downward to initiate spheroid transfer (**Fig. 6**).

**Note:** Ensure that each well of the ULA 384-well plate contains 35 - 40  $\mu$ L of cell culture medium for spheroid transfer. Excessive medium volume in the ULA 384-well plate may cause overflow after sandwiching the pillar plate onto the ULA 384-well plate. Avoid wetting the bottom of the pillars during this process, as it may lead to cross-talk or contamination.


**Note:** To position spheroids at the center of the pillars, ensure that there is no excess hydrogel on the pillars and gently insert the pillar plate into the ULA 384-well plate. Remove excess hydrogel from the pillars by horizontally sliding a 1 mL pipette tip across the pillar surfaces. In addition, carefully transport the inverted and sandwiched plates to a 5% CO<sub>2</sub> incubator without tilting, so as not to disturb vertical spheroid precipitation onto the pillars (**Fig. 8D**).



**Figure 6. (A)** Insertion of the 36PillarPlate into the ULA 384-well plate containing spheroids. **(B)** Inversion of the sandwiched plates for spheroid transfer.

- Incubate the inverted and sandwiched plates in a humidified 5% CO<sub>2</sub> incubator at 37°C for 4 hours to allow spheroid transfer from the ULA 384-well plate onto the pillar plate and attachment of the spheroids to the Matrigel.
- After incubation, carefully re-invert the sandwiched plates to return the 36PillarPlate to the top position, and then separate the pillar plate containing the transferred spheroids attached to the Matrigel.
- Immediately insert the pillar plate into the 384DeepWellPlate containing 70  $\mu$ L/well of prewarmed growth medium (**Fig. 7**).

**Note:** It is critical to prewarm the growth medium in the 384DeepWellPlate for at least 1 hour and gently tap the plate to dislodge any air bubbles in the wells before sandwiching with the pillar plate. Small micro-bubbles may appear at the edges of the pillars after sandwiching (**Fig. 8C**); however, these bubbles typically disappear within 1 - 2 days during routine medium changes.




**Figure 7. (A)** The 36PillarPlates sandwiched onto the 384DeepWellPlate. **(B)** Close-up image of the pillars of the 36PillarPlate inserted into the wells of the 384DeepWellPlate for cell culture.

- Inspect the pillar plate under a brightfield microscope to confirm successful spheroid transfer onto

the pillars (**Fig. 8**).

24. Culture the spheroids on the pillar plate in a humidified 5% CO<sub>2</sub> incubator at 37°C, replacing the culture medium every 1 - 2 days for culture using the 384DeepWellPlate.

**Note:** Cells on the pillar plate may also be cultured under dynamic conditions using a 36PerfusionPlate or petri dish combined with a digital rocker or low-speed rocker. However, low shear stress conditions (e.g., 5° tilt angle with 5-minute interval rocking) should be maintained for the first 4 days to prevent spheroid detachment from the pillar plate. Refer to the protocols titled “Dynamic Cell Culture in Perfusion Plate” or “Dynamic Cell Culture with PetriLid” for additional details.



**Figure 8.** (A) Stitched image of the entire 36PillarPlate with transferred spheroids. (B) Single pillar showing successful spheroid transfer. (C) Single pillar with micro-bubbles on the surface. (D) Single pillar with a spheroid located at the corner due to improper spheroid transfer.