

## **Beyond the Petri Dish: High-Throughput 3D Microbial Assays on the 384PillarPlate**

For more than a century, microbiology has largely been **two-dimensional (2D)**. Microbes are typically studied on flat agar plates or in liquid cultures. Yet in real environments, from dental plaque and chronic infections to industrial pipelines and water systems, microorganisms rarely exist in such simplified conditions. Instead, they organize into **complex three-dimensional (3D) communities known as biofilms**.

These biofilms behave very differently from free-floating microbial cells. They create protective extracellular matrices, develop nutrient and oxygen gradients, and often harbor **persister cells** that are far more resistant to antimicrobial treatments.

As microbiology research evolves, there is growing recognition that **3D microbial culture systems** are needed to better replicate real-world microbial behavior. One approach enabling this transition is the **384PillarPlate platform**, which allows researchers to rapidly generate hundreds of miniature 3D microbial cultures in a scalable and reproducible format.

### **Why the Third Dimension (3D) Matters**

Traditional antimicrobial assays expose microbes directly to drugs or preservatives in liquid culture or on flat agar surfaces. While these approaches are well established, they often fail to capture the structural and physiological complexity of microbial biofilms. Encapsulating microorganisms in a **three-dimensional (3D) matrix, such as agarose**, introduces several important features that more closely mimic natural microbial environments.

#### **Diffusion barriers**

In biofilms, antimicrobial compounds must penetrate through a protective extracellular matrix before reaching microbial cells. Similarly, when microbes are encapsulated in agarose microstructures, antimicrobial agents must **diffuse through the matrix**, creating more realistic exposure conditions.

#### **Nutrient and oxygen gradients**

Three-dimensional (3D) cultures naturally generate gradients of nutrients and oxygen. Cells located near the surface experience different metabolic conditions compared with those deeper within the matrix. These gradients often lead to the formation of **metabolically distinct subpopulations**, including persister cells that contribute to antimicrobial tolerance.

#### **Mechanical stability**

Surface-attached biofilms grown on plates can easily detach during washing or handling. In contrast, **agarose-encapsulated microbes on the pillar plate remain physically confined within the matrix**, maintaining structural integrity throughout experimental workflows.

Together, these properties make 3D microbial culture systems valuable tools for studying microbial physiology, antimicrobial resistance, and biofilm-associated behaviors.

## From Microbes to 384 Microcultures in Minutes

A major challenge in adopting 3D microbial culture systems has traditionally been **scalability**. Many biofilm models require specialized devices or lengthy preparation steps, limiting their use in high-throughput experiments.

The 384PillarPlate platform addresses this challenge through a simple workflow that enables rapid generation of hundreds of microbial microcultures.

### Step 1: Microbial encapsulation

Microorganisms are mixed with low-melting agarose to create a suspension that supports three-dimensional growth.

### Step 2: Stamping onto the pillar plate

The agarose-microbe mixture is transferred onto the 384PillarPlate through a straightforward stamping process. This step simultaneously loads microbes onto **all 384 pillars**, creating uniform micro-domes of encapsulated microbial cultures.

### Step 3: Growth and treatment

The loaded pillar plate is placed onto a complementary deep well plate containing growth media, antimicrobial agents, or other test compounds. Microbial viability and metabolic activity can then be measured using compatible readouts such as ATP bioluminescence assays.

Within minutes, researchers can generate **384 miniature microbial cultures on a single plate**, enabling rapid experimental throughput while maintaining controlled culture conditions.

## The Advantages of the 384PillarPlate

By combining 3D microbial encapsulation with a high-density pillar format, the platform enables scalable microbial assays with several practical advantages.

Feature	Advantage
<b>High throughput</b>	Up to 384 parallel microbial cultures in a single experiment
<b>Reproducibility</b>	Uniform microculture volumes and consistent growth conditions
<b>Speed</b>	Rapid microbial loading through a simple stamping workflow
<b>Efficiency</b>	Reduced reagent consumption compared with conventional assays
<b>Compatibility</b>	Works with standard plate readers, imaging systems, and automation tools

These characteristics make the platform particularly well suited for **screening-based microbiology workflows** where large experimental datasets are required.



## Applications of 3D Microbial Culture on the 384PillarPlate

While our recent work demonstrated the use of this system for **high-throughput antimicrobial efficacy testing**, the underlying methodology has much broader potential.

### Biofilm inhibitor discovery

Biofilms are notoriously difficult to treat due to their protective matrices and heterogeneous microbial populations. The 384PillarPlate enables rapid screening of compounds that disrupt biofilm structure or improve antimicrobial penetration.

### Personal care and household product testing

Preservatives and antimicrobial ingredients used in detergents, cosmetics, and personal care products must remain effective under realistic conditions. Miniaturized 3D microbial assays enable **rapid formulation screening while reducing reagent usage and experimental time**.

### Oral care research

Dental plaque represents one of the most common biofilms encountered in daily life. By tuning the density and growth conditions of encapsulated microbes, the platform can support screening of mouthwash formulations, toothpaste ingredients, and other oral care products.

### Industrial biofilm control

Biofilms are a persistent challenge in industrial environments such as water systems, food processing equipment, and pipelines. The 384PillarPlate platform enables high-throughput testing of disinfectants, coatings, and cleaning agents designed to control microbial growth.

### Antimicrobial resistance research

Three-dimensional microbial cultures naturally produce gradients that promote the formation of persister cells. This makes the system valuable for studying antimicrobial tolerance mechanisms and evaluating strategies to overcome biofilm-associated resistance.

### Microbiome and microbial ecology studies

Miniaturized microbial cultures allow researchers to evaluate how compounds influence microbial growth, metabolism, and community stability, providing opportunities to explore microbiome-friendly ingredients and microbial ecosystem dynamics.

## Bridging Laboratory Simplicity and Biological Complexity

Our recent publication in the *World Journal of Microbiology and Biotechnology* demonstrated the feasibility of using the 384PillarPlate platform for high-throughput antimicrobial efficacy testing using ATP-based viability readouts. By benchmarking the system against traditional colony-counting methods, we showed that miniaturized 3D microbial assays can provide **rapid and scalable alternatives to conventional antimicrobial testing workflows**. More broadly, the ability to generate hundreds of 3D microbial cultures in minutes opens new opportunities for microbial research across multiple industries.



## The Future of Microbial Screening

As microbiology continues to evolve, researchers increasingly require experimental systems that balance **biological realism with experimental scalability**. By combining **agarose-based 3D microbial culture with high-density pillar plate formats**, the 384PillarPlate platform provides a practical bridge between these two needs. The petri dish has served microbiology well for more than a century. However, the next generation of microbial research may well be **three-dimensional (3D), miniaturized, and high-throughput**.