

Selection guidance for pillar/perfusion plates and accessories

1. For static culture of 3D cells/organoids in a single condition

- 36PillarPlate (or 144PillarPlate)
- 36PetriLid (or 144PetriLid)
- LoadingPlate

2. For static culture of 3D cells/organoids in multiple conditions

- 36PillarPlate (or 144PillarPlate or 384PillarPlate)
- 384DeepWellPlate
- 36PetriLid (or 144PetriLid)
- LoadingPlate

3. For dynamic culture of 3D cells/organoids in a single condition

- 36PillarPlate (or 144PillarPlate)
- 36PetriLid (or 144PetriLid)
- LoadingPlate
- Equipment necessary: Any low-speed rocker in the lab

4. For dynamic culture of 3D cells/organoids in multiple conditions

- 36PillarPlate (or 144PillarPlate)
- 36PerfusionPlate (or 144PerfusionPlate)
- 36PetriLid (or 144PetriLid)
- LoadingPlate
- Equipment necessary: Digital rocker for bidirectional flow